

## Mycoplasma-ExS Efficacy Summary

| Activity Claim | Test Protocol | Contact Time | Test Organisms  |
|----------------|---------------|--------------|---|
| Bactericidal   | EN13697       | 5 minutes    | Pseudomonas aeruginosa, Escherichia coli, Enterococcus hirae, Staphylococcus aureus |

## Spore-EX Efficacy Summary

| Activity Claim | Test Protocol | Contact Time | Test Organisms  |
|----------------|---------------|--------------|---|
| Bactericidal   | EN1276        | 1 minute     | Escherichia coli, Staphylococcus aureus, Pseudomonas aeruginosa, Salmonella typhimurium, Enterobacter cloaca, Enterococcus hirae, Lactobacillus planetarium, Klebsiella pneumonia, Enterococcus faecalis, Proteus mirabilis, Streptococcus pyogenes, Salmonella choleraesuis, Listeria monocytogenes, |
|                | EN13727       | 1 minute     | MRSA, Klebsiella pneumonia, Staphylococcus aureus, Pseudomonas aeruginosa, Enterococcus hirae, Acinetobacter baumannii  |
|                | EN13697       | 1 minute     | Pseudomonas aeruginosa, Staphylococcus aureus, Escherichia coli, Enterococcus hirae   |
|                | EN1500*       | 30 seconds   | Escherichia coli  |
| Fungicidal     | EN1650        | 1 minute     | Candida albicans, Aspergillus niger, Tricophyton mentagrophytes   |
|                | EN13624       | 1 minute     | Aspergillus niger, Candida albicans   |
| Tuberculocidal | EN14348       | 1 minute     | Mycobacterium terrae, Mycobacterium avium   |
| Virucidal      | EN14476       | 1 minute     | HIV, Human Corona virus, Human Influenza virus, Herpes simplex virus, Swine Influenza A virus   |
|                | DVV/RKI**     | 1 minute     | Vaccinia virus, BVDV (Hepatitis C surrogate)  |
| Sporicidal     | EN13704       | 1 minute     | Clostridium difficile spores, Bacillus subtilis spores, Clostridium sporogenes spores, Bacillus cereus spores   |

\* The EN1500 test protocol (Chemical disinfectants and antiseptics. Hygienic handrub) was carried out on the hand wipe and hand foam presentations.

\*\*According to the DVV & RKI testing criteria, Spore-EX can be classified as active against all enveloped viruses in a 2% concentration in 1 minute.

### PromoCell GmbH

Sickingenstr. 63/65  
69126 Heidelberg  
Germany

Email: [info@promokine.info](mailto:info@promokine.info)

[www.promokine.info](http://www.promokine.info)

### North America

Phone: 1 – 866 – 251 – 2860 (toll free)  
Fax: 1 – 866 – 827 – 9219 (toll free)

### Deutschland

Telefon: 0800 – 776 66 23 (gebührenfrei)  
Fax: 0800 – 100 83 06 (gebührenfrei)

### France

Téléphone: 0800 90 93 32 (ligne verte)  
Téléfax: 0800 90 27 36 (ligne verte)

### United Kingdom

Phone: 0800 – 96 03 33 (toll free)  
Fax: 0800 – 169 85 54 (toll free)

### Other Countries

Phone: +49 6221 – 649 34 0  
Fax: +49 6221 – 649 34 40